


# cobas b 101 Overview

One system to screen, diagnose and monitor diabetes, dyslipidemia, and metabolic disorders at the Point-of-Care



- Convenient Point-of-Care diagnostic testing for HbA1c and lipid profile on one device with small footprint.
- Precise and reliable results comparable to lab.
- Fast turnaround time for results facilitates clinical and administrative workflows.
- Small sample size of capillary blood, venous whole blood, or plasma\* ~1-2 drops.
  - 2 µl for HbA1c
  - 19 µl for lipid panel
- Bi-directional connectivity with cobas IT 1000 data manager or other systems via standard protocol POCT-1A.
- Compliant with Canadian laboratory accreditation requirements and National Glycohemoglobin Standardization Program guidelines.
- Ideal for pharmacies, physician offices, out-patient clinics and satellite labs.
- For more information, please Contact us via [www.roclicanada.com](http://www.roclicanada.com).